

Bridging the Gap: Hospital-to-Home Polysomnography with Nox A1s

The demand for high-fidelity sleep data is outstripping in-lab capacity, leading to prohibitive wait times.¹ **The Nox SAS solution using the Nox A1s PSG System** provides a revolutionary way to simplify PSG testing: a less complex, portable system that delivers high quality signals which, when manually scored, show high levels of agreement with Level 1 PSG scoring of sleep stages.²

Meeting the “Gold Standard” Requirements

The Nox SAS solution using Nox A1s is designed to meet the rigorous guidelines set by the American Academy of Sleep Medicine (AASM)⁴ and similar guidelines for Level 2 unattended studies:¹

Feature	AASM Level 2 Requirement	Included in Nox A1s Level 2 studies
Channel Count	≥7 channels (EEG, EOG, EMG, ECG, Airflow, Effort, SpO ₂)	Meets requirements
Sleep Staging	Required (EEG, EOG, EMG)	Full PSG-based sleep staging supported (EEG, EOG, EMG via Nox A1s system when configured for SAS studies)
Respiratory	Airflow, Effort, SaO ₂	Pressure sensors, RIP belts, Wireless Pulse Oximetry
Cardiac	ECG/Heart Rate	Single lead modified Lead II
Position	Body position monitoring	Integrated 3D acceleration sensor



Reaching More Patients Where They Are: Clinical Implications & Resource Optimization

Level 2 in-home PSG offers the best of both worlds. It expands access beyond the limits of in-lab capacity, better supports patients who are reluctant or unable to sleep in a lab, and reduces the time staff spend on setup and cleanup. The result is a more comfortable experience for patients and a more efficient workflow for staff.

- » **Wait-Time Reduction:** Transitioning patients to Level 2 home studies frees up in-lab beds for more complex titration or neuro-sleep cases.
- » **Efficiency:** The "free-to-move" wireless design reduces setup time for PSG hookup in lab or at home to 15-20 minutes.
- » **Data Integrity:** A research study of the Nox SAS found that successful and reliable recordings were collected from approximately 85% of the 960 participants on the first attempt, and 88.6% overall after two attempts. The study also demonstrated that the EEG and EOG signals could be acquired without significant artifact.⁵

Versatile Hardware for Diverse Clinical Needs

The Nox A1s is a single hardware solution for Level I and II testing.

- » **Traditional EEG:** Supports standard scalp electrode placement for full sleep architecture and arousals monitoring.
- » **Nox SAS (Simplified Polysomnography):** An innovative forehead-based setup designed for rapid hookup and high patient compliance, eliminating the need for scalp electrodes and chin EMG while still recording EEG/EOG and facial EMG.

Find out how to get the precision of the sleep lab – everywhere. Talk to Nox about the SAS for full PSG studies in both lab and home settings.

Reference Verification: * AASM Manual V3 is the guiding clinical framework.⁴ Where applicable, local national guidelines (e.g., ERS, ESRS, AASM) should also be consulted.

References

1. Braun, M., Stockhoff, M., Tijssen, M. et al. A Systematic Review on the Technical Feasibility of Home-Polysomnography for Diagnosis of Sleep Disorders in Adults. *Curr Sleep Medicine Rep* 10, 276–288 (2024). <https://doi.org/10.1007/s40675-024-00301-z>.
2. Lee YJ, Lee JY, Cho JH, Choi JH. Interrater reliability of sleep stage scoring: a meta-analysis. *Journal of Clinical Sleep Medicine*. 2022;18(1):193–202. doi:10.5664/jcsm.9538
3. Punjabi NM, Kaplan PW, Margolick J, Aurora RN. 0319 A Simplified Bipolar Frontal Montage for Recording and Staging Sleep. *Sleep*. 2018;41:A122–A122. doi:10.1093/sleep/zsy061.318 2
4. Troester, M. M., Quan, S. F., & Berry, R. B. (2023). The AASM manual for the scoring of sleep and associated events: Rules, terminology and technical specifications (Version 3) [Manual]. American Academy of Sleep Medicine.
5. Punjabi NM, Brown T, Aurora RN, Patel SR, Stosor V, Cho JH, Helgadóttir H, Ágústsson JS, D'Souza G, Margolick JB. Methods for home-based self-applied polysomnography: the Multicenter AIDS Cohort Study. *Sleep Adv*. 2022 Apr 29;3(1):zpac011. doi: 10.1093/sleepadvances/zpac011. PMID: 35601080; PMCID: PMC9119085.

Disclaimer

The Nox SAS head cable, Nox SAS body cable and Nox A1s are CE-marked medical devices intended for clinical use at the direction of a trained healthcare professional. For more information on specific Nox Medical devices as part of the service offering, please consult manufacturer's Instructions for Use, including indications, contraindications, warnings and precautions at [noxmedical.com/downloads](https://www.noxmedical.com/downloads). Availability of Nox SAS may vary by region or market. For details, please contact Nox Medical or visit www.noxmedical.com