

Product Catalog

LBL-0102-REV03

Issue Date: 11th of September 2019



Contact

Nox Medical USA
5000 Research Ct
30024 Suwanee, GA
United States

Tel: (844) 475 9979

Fax: (678) 669 2274

sales@noxmedical.com

www.noxmedical.com

Nox Medical is a leading provider in the field of sleep diagnostic solutions that build on decades of experience in providing industry standard solutions for sleep monitoring and diagnostics as well as bleeding edge technology. Our main goal is to bring to market a new generation of products to support the growth of sleep medicine.


Nox Recording Devices have been engineered and designed to record and analyze sleep disorders in children as well as adults. With the monitors we have managed to remove as much sensor leads as possible since the patient comfort is the key element for good sleep recording.

Combining **Nox Recording Devices** with the powerful Noxturnal software unleashes the full potential of the Nox system. Offering automatic analysis, scoring and reporting, Noxturnal greatly simplifies the task of preparing a sleep study for the correct interpretation by the physicians.

Nox Medical will always aim to pro-actively satisfy its customers' needs. The customer service representatives are trained to handle customer requests promptly and professionally. Personalized customer service, along with prompt delivery, are major priorities in keeping commitments to customers.

NOX SYSTEMS	2
Nox T3® System with WristOx2, US	2
SOFTWARE	3
Licenses	3
Noxturnal EBM PG File Format License	3
Noxturnal EDF File Format License	3
SYSTEM EXTENSIONS	3
WristOx2™ Kit, Model 3150	3
Nox W7 Link Kit, S	4
Nox W7 Link Kit, R	4
Nox W7 Link Kit, A	5
CONSUMABLES	5
Nox RIP Belts	5
Nox RIP Belts Disposable, Pediatric	5
Nox RIP Belts Disposable, Small	5
Nox RIP Belts Disposable, Medium	5
Nox RIP Belts Disposable, Large	5
Nox RIP Belts Disposable, Extra Large	5
Cannulas	6
Oximeter Sensors	7
Leads and Electrodes	9
Nox Sensor Kits	10
Cleaning	11
SYSTEM COMPONENTS	11
Nox system components	11
Cross – Reference for Vyair customers	16

NOX SYSTEMS

REF NO	UNIT PR PACK	DESCRIPTION
51 11 10	1	Nox T3® System with WristOx2, US
		<p>The Nox T3 System is a high-end respiratory sleep diagnostic system. It comes with the Noxturnal analyses software. The Nox T3 is designed for both adult and pediatric patients.</p> 
THE NOX T3 SYSTEM INCLUDES		
56 11 10		Nox T3 Recorder
56 20 10		Nox Abdomen Cable
95 10 32		Nox RIP belts, Disposable M - 2 sets
55 20 10		Nox Cannula with Filter - 2pieces
56 20 11		Nox USB Cable
91 00 05		Nox T3 Carry Case
56 90 14		Nox Lid key
54 10 10		Nonin 3150 WristOx2™ Pack, includes
55 30 20		WristOx2 Soft Sensor, medium size
56 40 42		SPO2 3150 Wrist Band
91 00 08		CD-ROM of the Operators Manual
53 90 10		Noxturnal CD
91 10 01		Hookup Card
91 00 02		AA Battery - 2 units
91 00 01		AAA Battery - 2 units

WARRANTY COVERED FOR 2 YEARS (APPLIES ONLY TO THE MAIN HARDWARE UNIT)

SOFTWARE

LICENSES


REF NO	UNIT PR PACK	DESCRIPTION
53 61 10	1	Noxturnal EBM PG File Format License
<p>This license enables Noxturnal to read and review polygraph recordings in the Embla (EBM) file format.</p>		

REF	UNIT PR PACK	DESCRIPTION
53 61 11	1	Noxturnal EDF File Format License
<p>This license enables Noxturnal to read and review recordings in the EDF file format.</p>		


REF NO	UNIT PR PACK	DESCRIPTION
53 61 12	1	Noxturnal EBM PSG File Format License
<p>This license enables Noxturnal to read and review PSG recordings in the Embla (EBM) file format.</p>		

REF NO	UNIT PR PACK	DESCRIPTION
53 61 13	1	Noxturnal SomnoStar File Format License
<p>This license enables Noxturnal to read and review PSG recordings in the SomnoStar file format.</p>		


SYSTEM EXTENSIONS

REF NO	UNIT PR PACK	DESCRIPTION
54 10 10	1	WristOx2™ Kit, Model 3150
<p>Engineered with Nonin Medical's proven PureSAT SPO2 technology and updated with additional performance and reporting features, the WristOx2™, Model 3150 is the most advanced wrist-worn pulse oximeter available. The Nonin 3150 is an integrated part of the Nox sleep recording devices and provides Pulse Rate and Blood Oxygen Saturation signals.</p>		
INCLUDES		
55 30 20	WristOx2 Soft Sensor, medium size	
56 40 42	SPO2 3150 Wrist Band	
91 00 08	CD-ROM of the Operators Manual	


WARRANTY COVERED FOR 2 YEARS

REF NO	UNIT PR PACK	DESCRIPTION
54 40 10	1	Nox W7 Link Kit, S
		<p>The Nox W7 Link Kit - S uses Bluetooth technology to allow data transfer between the SenTec Digital Monitor and the Nox A1/T3 recorder for ambulatory studies.</p> 
		INCLUDES
		Nox W7 Link
		Gender Changer
		Nox W7 Power Adapter Kit
		with EU/UK/US Plug Adapters
56 20 11		Nox Mini USB Cable
		Leaflet

WARRANTY COVERED FOR 1 YEAR

REF NO	UNIT PR PACK	DESCRIPTION
54 40 11	1	Nox W7 Link Kit, R
		<p>The Nox W7 Link Kit - R uses Bluetooth technology to allow data transfer between supported auxiliary devices, such as the ResMed S9, Nonin RespSense Capnograph and selected Radiometer devices and the Nox A1/T3 recorder for ambulatory studies.</p> 
		INCLUDES
		Nox W7 Link
		Gender Changer
		Nox W7 Power Adapter Kit
		with EU/UK/US Plug Adapters
56 20 11		Nox Mini USB Cable
		Leaflet


WARRANTY COVERED FOR 1 YEAR

REF NO	UNIT PR PACK	DESCRIPTION
54 40 12	1	Nox W7 Link Kit, A
<p>The Nox W7 Link Kit - A uses Bluetooth technology to allow data transfer between supported auxiliary devices, such as the ResMed AirSense devices and the Nox A1/T3 for ambulatory studies.</p>		
INCLUDES		
	Nox W7 Link	
	Gender Changer	
	Nox W7 Power Adapter Kit	
	with EU/UK/US Plug Adapters	
56 20 11	Nox Mini USB Cable	
	Leaflet	

WARRANTY COVERED FOR 1 YEAR

CONSUMABLES

NOX RIP BELTS

REF NO	SETS PR PACK	DESCRIPTION
55 10 10	20	Nox RIP Belts Disposable, Pediatric
55 10 20	20	Nox RIP Belts Disposable, Small
55 10 30	20	Nox RIP Belts Disposable, Medium
55 10 40	20	Nox RIP Belts Disposable, Large
55 10 50	14	Nox RIP Belts Disposable, Extra Large
<p>Specifications</p> <ul style="list-style-type: none"> » Nox RIP Belts are very accurate in determining the effort of breathing » To be used with Nox Recording Devices » High quality signals » Super stretchable elastics » Comfortable for patient and clinician » Easy to adjust » Single use/disposable 		
		

NOT WARRANTY COVERED

CANNULAS

REF NO	UNIT OR PACK	DESCRIPTION
55 20 10	40	Nox Cannula with Filter, Adult

Specifications

- » The filter fits directly into the pressure connector on Nox recording devices
- » Total length is 90cm/35in
- » Air flow pressure
- » Single use/disposable



NOT WARRANTY COVERED

REF NO	UNIT OR PACK	DESCRIPTION
55 20 20	50	Nox Cannula with Luer-lock, Adult

Specifications

- » Male Luer Connector
- » Total length is 90cm/35in
- » Air flow pressure
- » Single use/disposable



NOT WARRANTY COVERED


REF NO	UNIT PR PACK	DESCRIPTION
55 21 51	25	Cannula with Luer-lock, Pediatric


Specifications

- » Male Luer Connector
- » Total length is 210cm/84in
- » Gas sampling line
- » Air flow pressure
- » Single use/disposable




NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
55 21 10	50	Nox Filter Tube Connector
<p>Connector between Nox Recording device and Nox Cannula with Luer-lock or Mask tubing.</p>		
NOT WARRANTY COVERED		

REF NO	UNIT PR PACK	DESCRIPTION
55 23 20	50	Mask tubing 183cm (72in) Female x Male
<p>The tube connects between the Nox Filter Tube Connector (Ref 55 21 10) and a CPAP Mask</p>		
NOT WARRANTY COVERED		

NOT WARRANTY COVERED

OXIMETER SENSORS

REF NO	SETS PR PACK	DESCRIPTION
55 30 10	1	WristOx2 Soft Sensor, Small
55 30 20	1	WristOx2 Soft Sensor, Medium
55 30 30	1	WristOx2 Soft Sensor, Large
<p>Specifications</p> <ul style="list-style-type: none"> » Preferred application: Fingers and toes » Reusable 		

WARRANTY COVERED FOR 2 YEAR

REF NO	UNIT OR PACK	DESCRIPTION
55 31 10	1	Flex Sensor with flexiwraps, Neonatal
<p>Specifications</p> <ul style="list-style-type: none"> » Patient Range < 2kg (4lbs) » Preferred application: Lateral aspect of foot » 1 meter/40 in cable » Reusable, made of durable silicone » Packed with 25 disposable FlexiWraps for neonatal <p>Note: Flex Sensors 3150I Sensor Adapter Cable is needed, Ref No 55 39 10.</p>		



WARRANTY COVERED FOR 3 MONTHS

REF NO	UNIT PR PACK	DESCRIPTION
55 31 20	1	Flex Sensor with flexiwraps, Infant
<p>Specifications</p> <ul style="list-style-type: none"> » Patient Range 2-20kg (4-44lbs) » Preferred application: Toe, thumb, medial aspect of hand » 1 meter/40 in cable » Reusable, made of durable silicone » Packed with 25 disposable FlexiWraps for infants <p>Note: Flex Sensors 3150I Sensor Adapter Cable is needed, Ref No 55 39 10.</p>		



WARRANTY COVERED FOR 3 MONTHS

REF NO	UNIT PR PACK	DESCRIPTION
55 31 30	1	WristOx2 Flex Sensor with flexiwraps, Adult
<p>Specifications</p> <ul style="list-style-type: none"> » Patient Range > 20kg (> 44lbs) » Preferred application: Index, middle or ring fingers » Reusable, made of durable silicone » 30cm/12 in cable » Ideal for numerous applications, including extended monitoring and sleep studies » Packed with 25 disposable FlexiWraps for adults <p>Note: Flex Sensors 3150I Sensor Adapter Cable is needed, Ref No 55 39 10.</p>		




WARRANTY COVERED FOR 3 MONTHS

REF NO	UNIT PR PACK	DESCRIPTION
55 31 11	25	FlexiWraps, Neonatal
55 31 21	25	FlexiWraps, Infant
55 31 31	25	FlexiWraps, Adult

Specifications


- » FlexiWraps to use with Nonin's Flex Sensors
- » Single use/disposable



NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
55 39 10	1	WristOx2 3150 Sensor Adapter Cable

All Nonin's sensors that are not WristOx2 sensors can be used with WristOx2TM, Model 3150 with this adapter cable.



WARRANTY COVERED FOR 3 MONTHS

LEADS AND ELECTRODES

REF NO	UNIT PR PACK	DESCRIPTION
		Nox Snap-on Lead:
55 40 20	2	- 50cm (20in), White
55 40 21	2	- 30cm (12in), Beige-White
55 40 22	1	- 100cm (40in), Green
55 40 23	1	- 50cm (20in), Beige-Green
55 40 24	2	- 150cm (60in), Grey
55 40 25	2	- 100cm (40in), Beige-Grey
55 40 26	2	- 150cm (60in), Black
55 40 27	2	- 100cm (40in), Beige-Black
55 40 28	2	- 100cm (40in), Orange

Specifications

- » Electrode cable with Female 1.5mm Touch proof safety connector and snap on head
- » For use with button snap type electrodes
- » White colored leads are in ideal lengths for Bruxism recordings
- » Green colored leads are ideal for the patient ground
- » Grey and black colored leads are ideal for LM recordings
- » Orange colored leads are ideal for ECG recordings
- » The leads with beige colored heads are ideal for pediatrics
- » Reusable



WARRANTY COVERED FOR 3 MONTHS

REF NO	UNIT PR PACK	DESCRIPTION
--------	--------------	-------------

55 42 09	25	Blue Sensor® N Electrode (Small)
----------	----	----------------------------------

Specifications

- » The small size makes it easy to use in the facial area and on pediatrics
- » 3.0 x 2.2cm (1.18 x 0.90in) snap on electrode with wet gel
- » Silver/silver chloride sensor
- » Adhesion ensures that the electrode stays in place
- » Single use/disposable



NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
--------	--------------	-------------

55 42 10	50	Blue Sensor® M Electrode
----------	----	--------------------------


Specifications

- » Ideal for EMG or ECG recordings
- » 3.4cm (1.25in) snap on electrode with wet gel
- » Silver/silver chloride sensor
- » Adhesion ensures that the electrode stays in place
- » Single use/disposable



NOT WARRANTY COVERED


CLEANING

REF NO	UNIT PR PACK	DESCRIPTION
55 90 10	1	Super Sani-Cloth PLUS Disinfection Wipes
<p>Specifications</p> <ul style="list-style-type: none"> » Disposable disinfection wipes for medical devices » Sani-Cloth is contraindicated for use as a terminal sterile/high level disinfectant on any surfaces 		


NOT WARRANTY COVERED

SYSTEM COMPONENTS


NOX SYSTEM COMPONENTS

REF NO	UNIT PR PACK	DESCRIPTION
56 20 10	1	Nox Abdomen Cable
<p>Connects the Abdomen RIP belt to the Nox recording devices.</p>		

WARRANTY COVERED FOR 6 MONTHS

REF NO	UNIT PR PACK	DESCRIPTION
56 20 51	1	Nox RIP Cable E, Abdomen
<p>Nox RIP cables E are designed to enable usage of Nox Respiratory Inductance Plethysmography (RIP) belts with 3rd party sleep recording devices. The cables along with Nox RIP Belts enable recording of respiratory effort signals. The Nox RIP cables have been validated for use with Embletta™ PDS and Gold portable sleep diagnostic systems.</p>		


WARRANTY COVERED FOR 6 MONTHS

REF NO	UNIT PR PACK	DESCRIPTION
56 20 52	1	Nox RIP Cable E, Thorax
		<p>Nox RIP cables E are designed to enable usage of Nox Respiratory Inductance Plethysmography (RIP) belts with 3rd party sleep recording devices. The cables along with Nox RIP Belts enable recording of respiratory effort signals. The Nox RIP cables have been validated for use with Embletta™ PDS and Gold portable sleep diagnostic systems.</p>
		

WARRANTY COVERED FOR 6 MONTHS

REF NO	UNIT PR PACK	DESCRIPTION
56 21 10	1	Nox RIP Cable Kit E
		A kit containing RIP cables E.
56 20 51	Nox RIP Cable E, Abdomen - 1 unit	
56 20 52	Nox RIP Cable E, Thorax - 1 unit	


WARRANTY COVERED FOR 6 MONTHS

REF NO	UNIT PR PACK	DESCRIPTION
56 90 10	1	Nox Service Kit
		<p>A kit containing replacement parts for the Nox recording devices.</p>
		
56 20 10	Nox Abdomen Cable - 1 unit	
56 90 11	Nox Battery Lid - Black	
56 90 13	Nox Clip Strap - 1 set / 2 clips	


WARRANTY INFORMATION GIVEN AS SINGLE ITEM

REF NO	UNIT PR PACK	DESCRIPTION
56 90 11	1	Nox Battery Lid (T3/A1)
Replacement Battery Lid for the Nox recording devices.		


NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 90 13	1	Nox Clip Strap, 2 pcs
Replacement Clips for the Nox recording devices.		

NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 90 18	1	Nonin 3150 Battery Lid
Replacement Battery Lid for the Nonin WristOx2, Model 3150		


NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 80 10	1	Nox T3 Carry Case
Replacement Carry Case for the Nox T3 System		


NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 90 14	1	Nox Lid key
Key to open the Nox device's battery Lid.		

NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 20 11	1	Nox Mini USB Cable
Used to connect the Nox recording devices to a computer. With standard USB on one end and mini USB on the other end.		


NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 40 42	1	WristOx2 3150 Wristband
<p>Specifications</p> <ul style="list-style-type: none"> » Wristband for the Nonin 3150 wrist pulse oximeter » Comfortable for both patient and clinician » Easy to adjust 		


NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 40 50	250	WristOx2 3150 Wristband, Disposable
<p>Specifications</p> <ul style="list-style-type: none"> » Single use wristband for the Nonin 3150 wrist pulse oximeter » Attaches the oximeter to patient for Chain of Custody » Easy to adjust 		

NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 30 91	4	MAHA Powerex Rechargeable 2700mAh Batteries
<p>Specifications</p> <ul style="list-style-type: none"> » High quality rechargeable batteries. » Nickel-metal hydride (NiMH) chemistry » 4 batteries to a pack. » 1.2V charge and 2700mAh capacity » Use with 56 30 90 Battery Charger MAHA 		

NOT WARRANTY COVERED

REF NO	UNIT PR PACK	DESCRIPTION
56 30 90	1	MAHA Powerex Battery Charger MH-C801D
<p>Specifications</p> <ul style="list-style-type: none"> » The battery charger Nox Medical recommends » Eight Cell 1-Hr PRO AA/AAA Charger » LCD Display » Worldwide heavy-duty power adapter 		

WARRANTY COVERED FOR 2 YEARS

CROSS - REFERENCE FOR VYAIRE CUSTOMERS

VYAIRE PART#	NOX PART#	DESCRIPTION
69446	56 90 10	Nox Service Kit
69447	55 21 10	Nox Filter Tube Connector
69590	55 10 20	Nox RIP Belts Disposable, Small
69591	55 10 30	Nox RIP Belts Disposable, Medium
69592	55 10 40	Nox RIP Belts Disposable, Large
69593	55 20 20	Nox Cannula with Luer-lock, Adult
69594	55 20 10	Nox Cannula with Filter, Adult
21816-001	55 10 10	Nox RIP Belts Disposable, Pediatric
21817-001	56 20 10	Nox Abdomen Cable
21818-001	56 80 10	Nox T3 Carry Case
22857-001	51 11 10	Nox T3 System with WristOx2, US
22988-001	56 11 10	Nox T3 Recorder
22989-001	56 40 42	WristOx2 3150 Wristband
22990-001	55 23 20	Mask tubing 183cm (72in) Female x Male
22991-001	55 21 51	Cannula with Luer-lock, Pediatric
22992-001	56 20 11	Nox Mini USB Cable
25635-001	56 90 13	Nox Clip Strap, 2 pcs
27486-001	56 90 11	Nox Battery Lid (T3/A1)
27641-001	55 10 50	Nox RIP Belts Disposable, Extra Large
27667-001	56 90 14	Nox Lid key

All information contained in this catalog is subject to change without notice.

Nox Medical is not responsible for any loss which results from any error in the product code or identification contained herein.

If you have any comments about our products and accessories, please help us to make our products even better by calling us or sending an e-mail: support@noxmedical.com.