



# Software and Firmware Update

## Release of Noxturnal 7.1

At Nox Medical, we are committed to continuous innovation and excellence. With the release of our latest Noxturnal 7.1 update, you will discover performance improvements and new advanced features, all crafted to enhance your experience with Nox Medical’s sleep diagnostic technology.

### Noxturnal 7.1 contains updates of the following

- Noxturnal Software
- Noxturnal Language Pack
- Nox T3s Firmware
- Nox A1s Firmware
- Nox C1 Firmware

### Main updates include:

- New feature: Support for 10-20 referential and bipolar EEG recording through the integration with Lifelines R40 Neuro amplifier.\*
- Improvements to 4% AHI reporting
- New feature: Ability to start a recording by attaching the Nox RIP belt.
- New feature: Successful Study Indicator.
- New feature: Ability to wake Nox A1s via Bluetooth.
- Important stability improvements and enhancements

**DISCLAIMER:** To take advantage of these new features, update to Noxturnal 7.1 and ensure your Nox A1s, T3s, and C1 are updated to the latest firmware. See the release bulletin for a compatibility guide. \*Requires Nox A1s and C1.

[Download Noxturnal 7.1](https://noxmedical.com/downloads-usa/)

<https://noxmedical.com/downloads-usa/>

[Download Release Bulletin](https://noxmedical.com/release-bulletin-usa/)

<https://noxmedical.com/release-bulletin-usa/>

### Software Download

To download the Noxturnal 7.1 software and Language Pack update, please visit our website.

### Release Bulletin

Within our Release Bulletin Hub, you can access a comprehensive document detailing the advanced features and enhancements included in the Noxturnal 7.1 software and firmware updates.