

Noxturnal 7.1.1

Release Bulletin



Version: 1.0 – LBL-0540

Release Date: 18th of December

Copyright© 2024

Nox Medical - All rights reserved

Nox T3s[®], Nox A1s[®], Nox C1 and Noxturnal[®] are manufactured by:

Nox Medical ehf

Katrinartuni 2

IS - 105 Reykjavík

Iceland

Website: www.noxmedical.com

nox medical

Copyright Notice

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language or computer language, in any form, or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without the prior written authorization from Nox Medical.

Overview

This Release Bulletin announces the release of the following Noxturnal software and manual.

- Noxturnal – **7.1.1 - 39309**
- Noxturnal Manual v. 4.4

The main aim of this release is to introduce Noxturnal 7.1.1 which focuses on addressing key feedback from the market since the latest Noxturnal release. Main improvements are towards scenarios with multi-site setup and large datasets focused on:

- Improvements of indexing the recording library to prevent Noxturnal freezing
- Improve handling of moving a bulk of studies between folders
- Study archiving improvements to prevent loss of studies during archiving process
- Handling network instability to reduce likelihood of lost scoring.

Updates to the Noxturnal Manual includes instructions for improved cybersecurity of the Noxturnal software.

Downloading Software

The software can be downloaded from different download locations and is available on the Nox Support Site.

Noxturnal & Noxturnal Language Pack

www.noxmedical.com/download

[Noxturnal - Using Dropbox](#)

Upgrade Recommendation and workarounds for Nox Users

Noxturnal Upgrade Wizard

Noxturnal 7.1.1

Noxturnal 7.1.1 has a built-in upgrade wizard for all users upgrading from 5.X or 6.0 versions of Noxturnal. This upgrade wizard prompts the user with selections on if they want to upgrade with old settings as a part of the install procedure or if they want a fresh install.

Older Noxturnal Versions

If users are working on an older version of Noxturnal (3.X and 4.X) it is recommended to do a fresh install and setup of the software since so many things have changed in recent releases and updating on top will create more problems than it solves.