



Software Updates

Noxturnal 7.1.1 Update

The latest Noxturnal release is here! Noxturnal 7.1.1 brings important improvements based on user feedback, especially for scenarios with multi-site setup and large datasets.

Key Enhancements Include:

- **Improved indexing:** Prevents freezing when navigating the recording library.
- **Better bulk study management:** Smoother handling when moving large numbers of studies between folders.
- **Enhanced study archiving:** Reduces the risk of data loss during archiving.
- **Stronger network stability:** Minimizes the chance of lost scoring due to connection issues.

Additionally, **Noxturnal Manual v. 4.4** now includes updated cybersecurity guidelines to help you protect your data.

[Download Release Bulletin](#)

Release Bulletin

For a detailed overview of the new features and improvements in version 7.1.1, check out the comprehensive document available in the Release Bulletin Hub on our website.

<https://noxmedical.com/release-bulletin-hub>

[Download Noxturnal 7.1.1](#)

Downloading the Software

The software can be downloaded from different download locations and is available on the Nox Support Site.

[Download Noxturnal 7.1.1 using Dropbox](#)

<https://noxmedical.com/downloads-row/>

<https://www.dropbox.com/scl/fi/8fs62tnrdh6eqrkp1ydv/Noxturnal-7.1.1.39309>